


# FULL SCORE FOR PROFESSIONAL SURFACE DISINFECTION

## Incidin™ Pro

The aldehyde-free concentrate for surfaces

Because of its SMART-BALANCED formulation, Incidin™ Pro brings a new standard to surface cleaning and disinfection.

- ▲ Your Pro for efficacy
- ▲ Your Pro for speed
- ▲ Your Pro for compatibility



# Incidin™ Pro

Full score for Professional surface disinfection

## Why Use Incidin Pro?

### BROAD SPECTRUM EFFICACY WITH SHORT CONTACT TIMES

Incidin Pro's antimicrobial efficacy is complimented by its short contact times and excellent material compatibility.

Effective against bacteria, yeast and various viruses (including norovirus), Incidin Pro is suitable for high risk areas and a wide range of surfaces, from general wards and public areas to operating rooms and intensive care units.

### EXCELLENT CLEANING EFFICACY AND MATERIAL COMPATIBILITY

It has been extensively tested on all surface types commonly found in the healthcare settings - from metals, linoleum and PVCs to ceramics, rubbers and plastics - with excellent results.

### EN STANDARDS COMPLIANT

Tests conducted to healthcare relevant EN standards confirmed the broad efficacy spectrum of Incidin Pro, including Norovirus, according to latest DVV and EN standard. See the table above for detailed test results.

### ALDEHYDE-FREE FORMULATION WITHOUT PERFUME AND DYE

The non-volatile and fragrance free ingredient system ensures a convenient and safe product use.

## MICROBIOLOGICAL AND VIROLOGICAL TEST RESULTS

EFFICACY	CONCENTRATION		CONTACT TIME
	%	g/L	
Disinfection of surfaces acc. to VAH with mechanical action (bactericidal, yeasticidal) (clean & dirty conditions)	0.5	5	30 min
	1.0	10	15 min
Disinfection of surfaces acc. to DGHM (German Society for Hygiene and Microbiology) method 2001 with mechanical action (tuberculocidal and mycobactericidal) (clean & dirty conditions)	2.0	20	2 h
	4.0	40	1 h
EN 16615 (bactericidal/yeasticidal) (clean & dirty conditions)	0.5	5	30 min
EN 13727 (bactericidal) (clean & dirty conditions)	0.25	2.5	5 min
EN 13624 (yeasticidal) (clean & dirty conditions)	0.25	2.5	5 min
EN 14348 (tuberculocidal, mycobactericidal) (clean & dirty conditions)	1.5	15	1 h
Effective against enveloped viruses (incl. HIV, HBV, HCV) acc. to RKI recommendation 01/2004*	0.25	2.5	30 min
	0.5	5	15 min
Adeno*	2.0	20	1 h
Papova*	1.5	15	15 min
Rota*	1	10	1 h
Rota*	1.5	15	5 min
Noro*	1.5	15	30 min
EN 14476 Noro (MNV)	1	10	1 h
	1.5	15	30 min

\* DVV 2015

### EUROPEAN HEADQUARTERS: ECOLAB EUROPE GMBH

Richtstrasse 7  
8304 Wallisellen  
Switzerland  
+41 (0) 44 877 2001  
www.ecolab.com

### ECOLAB DEUTSCHLAND GmbH

Ecolab-Allee 1  
40789 Monheim am Rhein  
Germany  
+49 (0) 2173 599 0  
www.de.ecolab.eu



# FAST. EFFECTIVE. ADVANCED.

## APPLICATION AREA

Incidin Pro is a liquid concentrate suitable for cleaning and disinfection of medical devices and all kind of surfaces in healthcare settings.

## INSTRUCTIONS FOR USE

Dilute concentrate with water (max 30°C) to the desired concentration. Wet entire surface with a sufficient amount of solution and leave to dry.

To achieve e.g. a 0.5% v/v solution prepare a 1:200 dilution with suitable water e.g. 5ml of Incidin Pro is added to 995ml of drinking water quality, giving a total volume of 1 litre.

## INGREDIENTS

100g contains the following active ingredients: 10.0g 2-Phenoxyethanol, 8.0g N,N-bis-(3-Aminopropyl) dodecylamine, 7.5g Benzalkoniumchloride.

## PACKAGING AND ORDERING INFORMATION

PRODUCT	OUTER SIZE	ORDER CODE
Incidin Pro sachet	400 x 20ml	XXXXXXX
Incidin Pro Bottle	3 x 2L	XXXXXXX
Incidin Pro Canister	6L	XXXXXXX
Incidin Pro Vat	30L	XXXXXXX
Incidin Pro Vat	200L	XXXXXXX

## STORAGE

Product should be stored in the original container below 25°C in a dry location out of direct sunlight.

USE BIOCIDES SAFELY. ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE.

