




Medline Industries Ltd
 No 5 Booths Park
 Chelford Road
 Knutsford – Cheshire
 WA16 8GS – England
 Phone: 0844 334 5237
 Fax: 0844 334 5238

www.medline.com/en/



Technical sheet

ORNEX Sterile Suction Tubing

Reference(s)

ORNEX56A, ORNEX66A, ORNEX76A, ORNEX510, ORNEX610, ORNEX710, ORNEX712

Description

Sterile suction tubing with 2 female ribbed connectors and 1 male connector pre-attached.

Code	Length (in meter)	Inner diameter (en mm)	Connection type*	Pack factor (each/case)
ORNEX56A	1,80M	5mm	F/F+M	50
ORNEX66A	1,80M	6mm	F/F+M	45
ORNEX76A	1,80M	7mm	F/F+M	40
ORNEX510	3,00M	5mm	F/F+M	30
ORNEX610	3,00M	6mm	F/F+M	25
ORNEX710	3,00M	7mm	F/F+M	20
ORNEX712	3,60M	7mm	F/F+M	20

***Connection type :**

F (female) connector, blue and inner ribbed, fits with male tapered connections for instrumentation/machines.

M (male) connector, white and transparent, is available in each pouch. It is pre-attached to female connector to receive female connector and can also be used to connect lengths of tubing together for longer lengths.

Main Materials

Tubulure: transparent PVC with DEHP
Blue connectors : Blue PVC // Male connector: transparent PVC
 Latex free

Pack factor / Packaging

Sterile double wrapped transparent individual peel pouch
 Pack factor : above

Sterilization

Sterilized by Ethylene Oxide

Manufacturer

Medline Industries, Inc - EN ISO 13485

Regulations

In compliance with Medical Device Directive 93/42/CEE – MDD Class: I sterile



© 2016 Medline Industries Inc. or one of its subsidiaries.
 All rights reserved.

TDS_ORNEX-
STERILE_EN02

01.01.INT.12