










# BD Vacutainer® Push Button Blood Collection Set



The **BD Vacutainer® Push Button** Blood Collection Set provides activation with a simple push of a button which ensures irreversible needle lock. The unique technology allows needle withdrawal from the vein and automatic retraction into the device, giving immediate protection against needlestick injury. The **BD Vacutainer® Push Button** Blood Collection Set allows for easy activation of the safety mechanism while attending to the patient/venipuncture site, and features:

- Safety activation confirmed by an audible click
- Blood flashback seen clearly through the translucent body
- Ease of use reducing training requirements

**BD Vacutainer®** EST tubes may be used as a discard tube, to remove the dead space from the tubing of blood collection sets. The BD EST tubes are plastic tubes and have no additive in them. See page 23 for details.

<b>Cat. No.</b>	<b>NHS Code</b>	<b>Gauge</b>	<b>Length</b>	<b>Shield Colour</b>	<b>Box/Case Qty</b>
<b>BD Vacutainer® Push Button Blood Collection Sets with Luer Adapter</b>					
367338	KFK324	21G	0.75", 7" tubing		50/200
367344	KFK363	21G	0.75", 12" tubing		50/200
367336	KFK323	23G	0.75", 7" tubing		50/200
367342	KFK366	23G	0.75", 12" tubing		50/200
367335		25G	0.75", 7" tubing		50/200
367341		25G	0.75", 12" tubing		50/200
<b>BD Vacutainer® Push Button Blood Collection Sets without Luer Adapter</b>					
367326		21G	0.75", 12" tubing		50/200
367324		23G	0.75", 12" tubing		50/200
367323		25G	0.75", 12" tubing		50/200