

Veterinary & Welfare

★ **Smart Oxygen** (Patent Pending)

new & improved!



A simple and effective way to provide oxygen or nebulising therapy, this is a cleverly designed PVC canopy which encloses an MDC cat basket and attaches to an oxygen source.

- ✓ **Instant oxygen or nebulising therapy**
- ✓ **Full patient visibility**
- ✓ **Includes plastic coated basket with base tray**

We are proud to announce exciting new developments to this 'must have' product. The Smart Oxygen is now available in two sizes; standard and (due to customer demand) a new large version. The PVC cover is now licensed as eco-friendly and we've added inlet and outlet scavenging ports with removable caps too.

Ideal for respiratory emergencies, pre and post GA oxygenation and nebulisation there is no need for restraint or masking and complete patient visibility may be maintained throughout the treatment.

No assembly required. Smart Oxygen includes plastic coated basket, base tray, bulkhead fittings and bag containing tubing, nebuliser, connector and 5 Venturi.



YOUR OPTIONS

- A236-360** Smart Oxygen - 46 cm x 29cm x 29cm
- A236-360LG** Smart Oxygen Large - 54cm x 38cm x 38cm



Assisting MDC with the testing of the new model was prominent critical care specialist Dr Marc Raffé DVM MS DACVAA DACVECC who said:

"The Smart Oxygen system is a practical and easy to use enclosure. When combined with the 'at a glance' ow rate guide , it provides life saving oxygen to all small animal patients."