

Blood collection sets

BD Vacutainer® Safety-Lok™ blood collection set




BD Vacutainer® Safety-Lok™ blood collection sets are sterile closed systems for venous blood collection. The safety mechanism is designed to help prevent needle stick injuries.

- Protection against needle injuries:
Following successful venepuncture, the integrated safety shield is pushed over the needle, surrounding it completely. It engages irreversibly with an audible click over the needle.
- Single hand activation possible:
The activation of the safety mechanism with a single hand allows greater attention to be paid to the patient and the venepuncture site.

- Versatile:
For taking blood samples and for short-term infusions of up to two hours.








BD Vacutainer® Safety-Lok™ blood collection sets without luer adapter

Cat No.	Size	Needle length	Tubing length	Luer adapter	Colour code	Unit of sales
367246	21G (0.8mm)	0.75" (19mm)	12" (305mm)	None		200 (4x50)
367247	23G (0.6mm)	0.75" (19mm)	12" (305mm)	None		200 (4x50)
368383	25G (0.5mm)	0.75" (19mm)	12" (305mm)	None		200 (4x50)



BD Vacutainer® Safety-Lok™ blood collection sets with luer adapter

Cat No.	Size	Needle length	Tubing length	Luer adapter	Colour code	Unit of sales
367282	21G (0.8mm)	0.75" (19mm)	7" (178mm)	With		200 (4x50)
367286	21G (0.8mm)	0.75" (19mm)	12" (305mm)	With		200 (4x50)
367284	23G (0.6mm)	0.75" (19mm)	7" (178mm)	With		200 (4x50)
367288	23G (0.6mm)	0.75" (19mm)	12" (305mm)	With		200 (4x50)
367295	25G (0.5mm)	0.75" (19mm)	7" (178mm)	With		200 (4x50)





Blood collection sets

BD Vacutainer® Safety-Lok™ blood collection set with pre-attached holder

With the pre-attached products, the holder is already fitted, so it is not necessary to manually assemble the needle and holder. The sterile closed system comes individually blister packed to minimise the risk of contamination of blood cultures. It is ideally suited for the taking of samples using the BD Bactec™ blood culture bottles.



BD Vacutainer® Safety-Lok™ blood collection sets with pre-attached holder

Cat No.	Size	Needle length	Tubing length	Colour code	Unit of sales
368654	21G (0.8mm)	0.75" (19mm)	7" (178mm)		200 (8x25)
368652	21G (0.8mm)	0.75" (19mm)	12" (305mm)		200 (8x25)
368655	23G (0.6mm)	0.75" (19mm)	7" (178mm)		200 (8x25)
368653	23G (0.6mm)	0.75" (19mm)	12" (305mm)		200 (8x25)