

Eurosonic 3D

The most balanced choice for ultrasonic washing



Ultrasonic washing tank with pre-set working temperature of 60°C and adjustable time up to 30 minutes. Eurosonic 3D offers maximum user-friendliness thanks to practical front buttons and luminous feedback on cycle status and correct operation. Its sturdy stainless steel tank (AISI 304L) and convenient side handles facilitate emptying, while the Sweep Mode technology ensures an effective and even cleaning.



Dimensions

27 x 24 x 19.5 cm (wxhxd)



Weight

3.7 Kg (with basket and cover)



Capacity

2.7 L



Supply voltage

220-230 V | 50/60 Hz

Power consumption

280 W

Accessories included



Anti-deformation non-drip tinted cover



Perforated stainless steel basket



0.6 L glass beaker with plastic lid



Cup holder 2 x 0.6 L



Fast

With Eurosonic Energy it is possible to start cleaning the instruments without the need to preheat the tub.



Simple to use

The cleaning cycle can be started with a simple gesture, in an intuitive way.



Intelligent

The indicator lights provide immediate information on the status of the cleaning cycle or on the presence of any errors or anomalies in its use.

