

# Safety blood collection sets

## BD Vacutainer® Safety-Lok™ blood collection set with pre-attached holder

With this safety blood collection set the holder is already pre-attached, so it is not necessary to manually assemble the needle and holder. This ready-for-use blood collection set is individually packaged in a sterile blister.





The sterile closed system is ideally suited for the taking of samples using the BD BACTEC™ blood culture bottles.

### Further information

Technical information is available on request.



## BD Vacutainer® Safety-Lok™ blood collection sets with pre-attached holder

Cat. no.	Size	Needle length	Length of tubing	Colour code	Box/case
368654	21G (0.8mm)	19mm	178mm		25/200
368652	21G (0.8mm)	19mm	305mm		25/200
368655	23G (0.6mm)	19mm	178mm		25/200
368653	23G (0.6mm)	19mm	305mm		25/200